www.zimed.com.tr

zimed®

**EKSTERNAL FİKSATÖRLER
KULLANIM KILAVUZU**
*INSTRUCTIONS FOR USE
FOR EXTERNAL FIXATORS*

TR

ENG

ES

IT

FR

RU

innovation to move

(TR) ÜRÜN TANIMI:

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ZİMED MEDİKAL

INSTRUCTONS FOR USE FOR EXTERNAL FIXATORS

Submitted to the Attention of Surgeons (EN)

Document No: IFU.05.01/00
14.06.2022

(EN) PRODUCT DEFINITION:

External fixators have wide-spread use in the treatment of fractures and complications, particularly of open and fragmented fractures, infected and non-infected nonunion and deformities.

External fixation systems are designed to be used with pins and Schanz screws in open fractures, intra- and extra-articular fractures for minimal fixation, fractures in poly-traumas, elongation of legs, osteotomy, arthrodesis and other indications concerning bones that repair is possible in adults and pediatric patients. External Fixators are offered both as sterile and non-sterile. Two different sterilization methods, gamma irradiation and ethylene oxide, are used for sterilization.

External fixation systems are into groups to be used in:

- Long bone fractures,
- Articular and peri-articular fractures,
- Pelvic fractures,
- In fractures of small bones of hands and feet.

Extremity Elbow Fixator: It is used in fractures of the elbow area for minimal invasive applications. It can be angulated and positioned around the elbow. It has standard measures. It is made from AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Segmental Fixator (Humerus): The product is locked onto the middle clamp fixator rail after proximal or distal osteotomy in humerus. Lengthening is carried out by moving each of the other proximal and distal clamps in their respective directions. There are 3 different length options as the short, medium and long. They are made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Dynamic Axial External Fixator (Humerus): Dynamic Axial External Fixator is used for stabilization in the humerus area, treatment of fractures, and treatment of multiple fractures in the upper extremity and femoral-tibial areas of pediatric patients. Multidimensional properties of the fixator allow simple angular adjustments. It is a mono-plan external fixator in ball-joint structure with connections with Schanz retaining clamps with intertwining telescopic body. There are 3 different length options as short, medium and long. It is made from AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

T-Clamp Humerus: It is used in proximal humeral fractures (in children and adolescents). It is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Wrist Fixator (Penning Type): The Wrist Fixator with dynamic angulated clamp can be used in the partial osteotomy for distal radius deformity, for acute repair by filling the space with graft or gradual repair of the space by callus tissue after hemicallotaxis osteotomy. It can be used in distal radius and/or ulna fractures, in instable intra-articular distal radius fractures (by fixing the fragments of the fracture with other implants). Lockable angulations can be carried out thanks to the ball-joint structure. Schanz screws can be placed in plans with different angles thanks to the clamps that angulation is possible. In distal radial fractures, it can be used with the application of other internal fixation implants and fixation of the joint. In can also be used in extra-articular fractures of the wrist, radial distal and shaft fractures and in ulnar fractures.

Hand Wrist Fixator has 3 different lengths as short, medium and long sizes. It is made of AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Finger Fixator (Joint): Finger Fixator is used for the fixation of fractures in finger joint areas. It has standard measures. It is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Dynamic Finger Fixator: Dynamic Finger Fixator is used for the fixation of fractures in finger bones. It has standard measures. It is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Lengthener Finger Fixator: Lengthener Finger Fixator is used for the fixation of fractures of finger bones and for lengthening the finger. There are 3 different length options as the short, medium and long. It is made from AISI 304, 6013 aluminum and TiAl4V Elv materials.

Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Segmental Fixator (Femur-Tibia): It is locked onto the middle clamp's fixator rail in very short femur or tibia for single- or double-focus lengthening after proximal and distal osteotomy. Lengthening is carried out by moving each of the other proximal and distal clamps in their respective directions. Osteotomy can be applied on the mid-shaft of the short femur and tibia for lengthening. In large proximal or distal bone losses and/or length losses of femur or tibia, proximal metaphyseal osteotomy can be performed to move the mid-segment to the distal by means of the clamp in the middle to ensure uniting. In the large, central injuries of femur or tibia with bone loss, proximal and distal metaphyseal osteotomy is performed. Simultaneously, proximal and distal transport will be carried out till the segments come together. In large peripheral injuries of femur and tibia with bone loss, two different osteotomy procedures will be performed on the long bone segment. Simultaneously, proximal movement of both mid-clamps will be carried out till the space is closed. There are 3 different length options as the short, medium and long. It is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Limb Lengthener Fixator: In segmental fractures of distal femur, in large unstable (irregular) tibial fractures, in fixation of distal tibial fractures with T-clamp for the proximal segment; in fixation of proximal tibial fractures with T-clamp in metaphysis, and in gradual angular correction of proximal tibia varus or valgus deformities. (The T-clamp is mounted in the proximal in anterior direction for osteotomy). In short tibia, for the gradual angular correction of proximal varus or valgus deformities, it can be used with the purpose of elongation on the deformity line. There are 3 different length options as the short, medium and long, and it is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Dynamic Axial External Fixator (Femur - Tibia): Dynamic Axial External Fixator has the capability of multi-directional and motile length adjustment for uses in articular stabilization between femoral and tibial areas, treatment of fractures, treatment of multiple fractures of lower extremities and use in distal tibial area. There are 3 different length options as the short, medium and long. It is made from AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

135-Degree Clamp: It is used in petrochanteric fractures of the femur. It is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

T-Clamp Femur: It is used in distal femoral fractures and valgus deformities and distal and proximal tibial fractures, varus deformities and osteotomy procedures. It is made from AISI 304 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Dynamic Axial External Fixator (Pelvic): Pelvic Dynamic Axial External Fixator is used for precise fixation of vertical fractures in pelvic area and correction of instability. There are 3 different length options as the short, medium and long. It is made from AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Pylon Fixator: Pylon fixator has been developed for easy and quick applications for fractures and deformities in the ankle joint area. It is a fixator specifically designed for tibial pylon fixator and ensures perfect fixation and allows motility of the joint in the early period. There are 3 different length options as the short, medium and long, and it is made from AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Hand-Wrist Fixator: It is used in the partial osteotomy of distal radial deformity by filling the space with graft followed by acute correction or hemicallotaxis osteotomy followed by gradual correction of the space with callus; in distal radial and/or ulnar fractures, and in unstable intra-articular distal radial fractures (through fixation of the fracture fragments with other implants). It has a standard length and it is made from AISI 304, 3Cr13 and 6013 aluminum materials. Aluminum materials are suitable for eloxal coating, and eloxal coating procedure is applied in this type of materials.

Schanz Screw (Cortical, with threads): It is used for fixation in places where bones are cortical structure. Schanz screws with cortical threads are made from materials manufactured according to ISO 5832-1 standards. Schanz Screws have 6 length options in total in 35, 40, 50, 60, 70, 80, 90, 100 mm lengths at Ø2 mm. There are 7 length options in total at Ø3 mm between 60 mm and 110 mm with sequential increments of 10 mm. There are 8 length options in total at Ø4 mm between 70 mm and 140 mm with sequential increments of 10 mm. There are 7 length options in total at Ø5 mm between 140 mm and 260 mm with sequential increments of 10 mm. There are 7 length options in total at Ø6 mm between 140 mm and 260 mm with sequential increments of 10 mm.

INDICATIONS:

External Fixator system is used for quick stabilization of open and/or instable fractures. It is ideal for tibia, femur and humerus, and the External Fixator system with multi-directional effects can be used in many fracture types, and in addition, in compression arthrodesis, extremity lengthening and bone segment transport cases.

Fixators are medical devices to immobilize fractures and they ensure stability in open or closed partial fractures. The external fixator system includes various components used to build structures that treat the indicators listed below. Different modular frame configurations are possible with external fixation systems. Individualized configurations must be designed for specific applications in each situation. Components of the entire external fixation system are always for single use. Each external fixation system is designed as a complete system and replacement with components of other systems is not allowed.

External Fixation Systems are made from aluminum and stainless steel. Circular frames and half-pin fixation techniques are used for the placement of the device. The device is designed at a central position on the rotational axis. Indications of external fixators must be considered in well-selected cases, where they are more advantageous as compared to conventional methods such as open reduction-internal fixation and casting.

Indications are divided into 3 groups including,

- A) Approved indications,
- B) Potential indications, and
- C) Controversial indications.

A) Approved indications:

- 1) Gustilo Anderson type 2 and type 3 open fractures (TYPE 2: Medium-level soft tissue injury larger than 1 cm; minimal periosteal separation; Type 3 is handled in 3 divisions: TYPE 3A: Severe tissue damage and abundant contamination, covering the wound with skin is possible; TYPE 3B: Severe tissue damage and abundant contamination, covering the wound with skin is not possible; TYPE 3C: open fracture with arterial injury requiring repair),
- 2) Fractures accompanied by burns,
- 3) Subsequent counter-lateral extremity flap (tissue harvested with preservation of blood supply), temporary immobilization in fractures that will require Free Tissue Flap or reconstruction,
- 4) Fractures that will require distraction later such as fractures with bone loss,
- 5) Extremity lengthening,
- 6) Arthrodesis procedure,
- 7) Pseudoarthrosis or infected fractures.

B) Potential indications:

- 1) Pelvic fractures and dislocations,
- 2) Open and infected pelvic pseudoarthrosis cases,
- 3) Reconstructive pelvic osteotomy (such as bladder exstrophy),
- 4) After complete removal of the tumoral tissue and application of auto graft or allograft,
- 5) After femoral osteotomy in children (in children, after internal fixation with plaque-screw, to remove the latter and exclude the 2nd operation),
- 6) To repair and restructure vascular and nerve injuries accompanying fractures,
- 7) Re-implantation of the amputated extremity,
- 8) Fixation of closed fractures in multiple extremities. Treatment with external fixators is the alternative treatment in patients with polytrauma. In these patients, fractures can be treated with traction, casting or open reduction. However, immobilization with these methods is difficult.
- 9) Correction of congenital joint contractures,
- 10) To support the internal fixation between major fragments in fragmented fractures,
- 11) Ligamentotaxis,

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlık Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

ZİMED MEDİKAL

INSTRUCTONS FOR USE FOR EXTERNAL FIXATORS

Submitted to the Attention of Surgeons (EN)

Document No: IFU.05.01/00
14.06.2022

12) To immobilize the fracture in patients with head trauma: Convulsions and continuous contractions occur in patients with severe head trauma. In such cases, fixation of the fractures with casting, traction or non-external fixation methods, and rigid fixation can fail in cases where convulsions are frequent and contractions are strong. With this reason, temporary external fixation can be performed to immobilize the fracture.

13) It can be used for temporary fixation in cases where the patient must be transferred with diagnostic or therapeutical purposes or because of other surgical conditions. External fixation will prevent the impairment of fracture alignment during transport.

14) It is very difficult to evaluate the knee ligaments in proximal tibial or distal femoral fractures. External fixator will allow the evaluation of knee ligaments if applied very close to the fracture.

C) Controversial indications:

Closed fractures: Application of external fixators in closed fractures and the successful conventional treatment are controversial. Even when meticulous attention is given in the basic principles of re-positioning, potential problems including infections in nail trace, delays in union and re-fracturing can occur. External fixation in the treatment of long bones is very beneficial. However, it must be used for patients, for whom reduction or immobility cannot be ensured with traditional methods.

Do not use this product rather than the indications stated here.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

CONTRA-INDICATIONS:

Active infections, manifest or occult infection or manifest local inflammation in or around the affected site, existing fracture fixation device implanted before, hypovolemia, hypothermia, coagulopathy, obesity, mental disorders preventing cooperation with the rehabilitation regime, advanced osteoporosis, poor bone quality or amount, disorders in the medullar canal, poor blood circulation, and previous infections will delay healing.

This product is contra-indicated in patients who are unwilling or incapable as regards in following the instructions for postoperative care, or who have mental or physiologic disorders. It is not recommended in patients with severe osteoporosis, who are HIV-positive or with immunosuppression, or poorly-controlled diabetes mellitus. Tests must be carried out before placement of the implant in patients with foreign body hypersensitivity or suspected sensitivity against materials.

WARNINGS:

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
3. Products of different manufacturers must not be used together.
4. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
5. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning, and must inform the patient.
6. These products are for single use only. Do not re-use the product.
7. Products must never be cleaned with ultrasonic washing detergents. Ultrasonic washing detergents can cause removal of the coating in external fixator products with eloxal coating.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:

Preoperative:

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection:

weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.

2. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.

3. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension, and can result in breaking of the product.

4. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.

5. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used.

6. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.

2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.

3. During the implantation, fractures can occur in the application area because of trauma or overloading, particularly in the presence of poor bone structure.

4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.

2. Abrasion of the components can result in the increase of loosening and damaging of the bone.

3. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending of cracking of components.

4. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.

5. Recommendations of surgeons in the postoperative period are of utmost importance.

6. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

SIDE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarcts can occur.
2. Damaging of blood vessels or hematomas can be seen.

3. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.

4. Unwanted lengthening or shortening of the extremity can occur.

5. Fatigue or breaks of products can happen as a result of trauma, exhausting activity, wrong alignment, defective placement of implants, period of service, loss of fixation, lack of union or excessive weight

6. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.

7. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.

8. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.

9. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT:

The period of the implant within the patient's body is minimum 6 weeks and can be extended to 1 year in difficult situations or more; afterwards, it must be removed based on the surgeon's decision.

STERILIZATION:

Implants Put in the Market as Sterile: External fixation products are put in the market as sterile or non-sterile. Sterile products will remain sterile until the package integrity is broken. There are two types of sterilization as sterilization with gamma irradiation or ethylene oxide. Sterilization method will be shown on the packaging. Minimum 25 kGy gamma irradiation is applied to all the implants sterilized with gamma irradiation.

Please inform the manufacturer about any products with impaired sterility although offered as sterile.

Implants Put in the Market as Non-Sterile: Products put in the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in.

All the intramedullary nail products and surgical tools used with them must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

Autoclave Sterilization: Non-sterile products put in the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

AUTOCLAVE STERILIZATION	
Parameters	Autoclave with pre-vacuum
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying Time	15 minutes

Note: Exposure time means the period to be applied starting from the moment that the autoclave temperature reaches the exposure temperature.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

ZİMED MEDİKAL

INSTRUCTONS FOR USE FOR EXTERNAL FIXATORS

Submitted to the Attention of Surgeons
(EN)

Document No: IFU.05.01/00
14.06.2022

The table below shows autoclave types for specific sterilization parameters.

PRODUCT TYPE	RECOMMENDED STERILIZATION METHOD
External Fixators	High Temperature Prevac Steam (Pre-Vacuum Autoclave) Flash Prevac Steam (Pre-Vacuum Autoclave): <i>If used, practice must be only according to ANSI/AAMI ST79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)</i>

Since the method of use will change for each device, please consult the written instructions of your manufacturer for your sterilizer and loading configuration.

SHELF LIFE:

Shelf life for **External Fixators** sterilized with gamma irradiation is 5 years; while the same for those sterilized using ethylene oxide is 3 years.

Please check the expiry date. Contact the manufacturer for products expired.

Non-sterile products do not have shelf lives.

DISPATCH OF THE PRODUCTS:

Intramedullary Nails presented as sterile or non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrank together with instructions for use and label.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity. The recommended storage conditions for sterile products are 15-35°C.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

(See: Description of Symbols)

SYMBOL DESCRIPTIONS	
	Logo of the Firm
	"1984" Notified Body Number
	Medical Device
	Lot Number
	Product Reference /Catalog Number
	UDI Number (Single Device Identification Number)
	Electronic Instructions for Use

	Manufacturer Details
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with irradiation
	Sterilized with ethylene oxide
	Used only with prescription
	Implants are for single use only
	Do not re-sterilize
	Do not use the damaged package
	Keep away from direct sunlight
	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution
	Non sterile

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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