
www.zimed.com.tr

zimed®

FEMORAL STEMLER
KULLANIM KILAVUZU
FEMORAL STEMS
USER GUIDE

innovation to move

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Attention of the Specialist (EN)

Document No: IFU.07.02/02
17.11.2025

IMPORTANT INFORMATION - Please read before use

Visit <https://zimed.com.tr/e-ifu> for the latest eIFU and access to different languages.

PRODUCT DESCRIPTION AND INDICATIONS:

(Do not use the products except for the indications written here)

DEVICE DESCRIPTION

Zimed Medical Femoral Stems consist of Cemented Femoral Stems, Cementless Femoral Stems, Cemented Calcar Replacement Stems, Cementless Calcar Replacement Stems, Plugs, and Distal Centralizers.

Femoral Stem provides support to the body by replacing the bone in cases of fractures or tumors in the femoral head and neck. It has an anatomical design.

Cementless Femoral Stem is manufactured from Ti6Al4V material in accordance with ISO 5832-3 standard and has no right-left distinction. For cementless applications, HA, Dual, and Porous coatings are applied. In cementless applications, the neck angle is 132°. The Femoral Stem is available in 7 sizes, ranging from 7.5 mm to 13.5 mm.

Cemented Femoral Stem is manufactured from CoCrMo according to ISO 5832-4 (ASTM F75) standard and from implant-grade steel according to ISO 5832-1 (ASTM F138) standard. It has no right-left distinction. In cemented applications, the neck angle is 130°. The Femoral Stem is available in various lengths, ranging from 7 mm to 12 mm.

Calcar Replacement Stem provides support to the body by replacing the bone in cases of fractures or tumors in the femoral head and neck. It has an anatomical design.

Cementless Calcar Replacement Stem is manufactured from CoCrMo material in accordance with ISO 5832-4 standard and has no right-left distinction. For cementless applications, HA, Dual, and Porous coatings are applied. In cementless applications, the neck angle is 135°.

Cemented Calcar Replacement Stem is manufactured from CoCrMo according to ISO 5832-4 (ASTM F75) standard and from implant-grade steel according to ISO 5832-1 (ASTM F138) standard. In cemented applications, the neck angle is 135°. It replaces the bone in cases of fractures and tumors in the femoral head region and provides support to the body.

Plug prevents cement leakage inside the bone and helps the stem anchor to the bone. It is manufactured from UHMWPE material in accordance with ISO 5834-2 standard and is available in 9, 10, 11, 12, 13, 14, 15, 16, 17, and 18 mm sizes.

Distal Centralizer centers the stem in the bone to correctly position the prosthesis. It is manufactured from UHMWPE material in accordance with ISO 5834-2 standard and is available in 8.5, 11.5, and 13.5 mm sizes. Suitable stems for each centralizer size are defined in the table below.

Distal Centralizer	Cemented Femoral Stem (Size)	Cemented Calcar Replacement Stem (Size)
8,5 mm	7 mm, 8 mm	9 mm, 10 mm, 11 mm
11,5 mm	9 mm, 10 mm	12 mm, 13 mm, 14 mm
13,5 mm	11 mm, 12 mm	15 mm, 16 mm, 17 mm

PROSTHESIS COATINGS

The dual coating (Ti + HA) and porous coating applied to the implant surface enhance implant stability by supporting osteointegration. The dual coating consists of titanium, which provides mechanical fixation, and hydroxyapatite, which promotes bone growth. The porous coating facilitates cellular infiltration, strengthening long-term biological fixation.

INTENDED USE

Zimed Medical Femoral Stems are used in the hip joint for conditions such as osteoarthritis, hip dysplasia, femoral head fractures, primary coxarthrosis (hip osteoarthritis), traumatic arthritis, slipped capital femoral epiphysis, fused hip, rheumatoid arthritis, congenital hip dislocation, avascular necrosis, and joint damage resulting from hip fractures.

INDICATIONS

- Osteoarthritis (Degenerative joint disease) (When the cartilage in the joint wears out due to aging)
- Rheumatoid arthritis (An autoimmune disease that causes chronic inflammation and damage in the joint)
- Post-traumatic arthritis (Secondary arthritis that develops after hip fractures or dislocations)
- Osteonecrosis of the femoral head (Avascular necrosis) (Bone death in the femoral head due to disrupted blood supply)
- Hip dysplasia (Congenital or developmental) (Abnormal development of the hip joint, leading to early cartilage damage with aging)
- Ankylosing spondylitis (A progressive inflammatory disease affecting the hip, leading to restricted movement and pain)
- Presence of benign or malignant tumors around the hip
- Inadequate response to conservative treatments (In cases where painkillers, physical therapy, injections, and similar treatments are ineffective)

Do not use the product except for the indications written here.

CONTRAINDICATIONS

- Conditions that prevent or tend to prevent adequate provision of implant support or prevent the use of an appropriately sized implant, e.g:
 - Limited blood supply;
 - Inadequate quantity and quality of bone support, e.g. osteoporosis or metabolic disorders that impair bone formation and osteomalacia;
 - Infections or other conditions leading to advanced bone resorption.
- Mental and nervous disorders that impair the patient's ability and willingness to limit activity
- Patients with active or suspected infection around the hip joint,
- Physical conditions and movements that place excessive stress on the implants. For example, Charkot joints, muscle deficiencies, multi-joint capability disorders, etc.
- Skeletal underdevelopment,
- Patients with any vascular insufficiency, muscular atrophy or severe neuromuscular disease,
- Patients with known moderate or severe renal impairment,
- Patients with suppressed immune systems due to diseases such as AIDS or high doses of corticosteroids,
- Physical ailments such as morbid obesity or restricted movement,
- Patients not following postoperative care instructions,
- Contraindications may be relative or absolute, and the patient should be carefully examined and alternative procedures performed if possible. For example, non-operative treatment, osteosynthesis, arthrodesis, femoral osteotomy, pelvic osteotomy, resection arthroplasty, semi-arthroplasty and others.
- High-risk conditions include: osteoporosis, metabolic disorders that prevent bone formation and osteomalacia.

PATIENT POPULATION

Designed for use in adult patients with a height of 145–190 cm and a body mass index (BMI) of 34 or below. General principles should be applied in patient and implant selection. Correct implant choice is very important. The appropriate type and size should be determined by considering anatomical and biomechanical factors such as the patient's age, activity level, weight, bone and muscle condition, and whether they have undergone previous surgeries.

TARGET USER

The intended user of Zimed Medical Femoral Stems is surgeons. These devices should not be used by unqualified personnel. It is very important that the surgeon and operating room

staff are fully familiar with the proper surgical technique for the use of the implant and related instruments.

RAW MATERIAL INFORMATION

Ti6Al4V ELI raw material is produced in accordance with ISO 5832-3 and ASTM F136. The raw material content is given in the table below.

Element	% Composition
Nitrogen (N), most	0.05
Carbon (C), most commonly	0.08
Hydrogen (H), most commonly	0.012
Iron (Fe), most commonly	0.25
Oxygen (O), most commonly	0.13
Aluminum (Al)	5.5-6.50
Vanadium (V)	3.5-4.5
Titanium (Ti)	90.48-88.48

CoCrMo raw material is produced in accordance with ISO 5832-4 and ASTM F75 standard. Raw material content is given in the table below.

Element	% Composition
Chromium (Cr)	27.0 – 30.0
Molybdenum (Mo)	5.0 – 7.0
Iron (Fe)	0.75 max.
Manganese (Mn)	1.0 max.
Silicon (Si)	1.0 max.
Carbon (C)	0.35 max.
Nickel (Ni)	0.50 max.
Nitrogen (N)	0.25 max.
Cobalt (Co)	63.7-58.71
Tungsten	0.20 max.
Phosphorus	0.020 max.
Sulfur	0.010 max.
Aluminum	0.10 max.
Titanium	0.10 max.
Boron	0.010 max.

Implant Steel (SS) raw material is produced in accordance with ISO 5832-1 and ASTM F138 standards. Raw material content is given in the table below.

Element	% Composition
Carbon (C), most commonly	0.030
Silicon (Si), most commonly	1.00
Manganese (Mn), most commonly	2.00
Phosphorus (P), most commonly	0.025
Sulfur (S), most	0.010
Nitrogen (N), most	0.10
Chromium (Cr), most commonly	17.0-19.0
Molybdenum (Mo), most commonly	2.25-3.0
Nickel (Ni), most commonly	13.0-15.0
Copper (Cu), most commonly	0.50
Iron (Fe), most commonly	35.92-60.67

In case of any negative situation with the product, please contact the manufacturer. Manufacturer information detail is given below.

PRODUCER: ZİMED Medikal San. ve Tic. Ltd. Şti.
Aydınlı Mah. 03070 Cad. No:4 Şehitkamil 27580 Gaziantep-Türkiye
Tel : 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:**www.zimed.com.tr **E-mail:** info@zimed.com.tr

UHMWPE raw material is produced in accordance with ISO 5834-2 standard. Raw material content is given in the table below.

Element	Maximum Amount (mg/kg)
Ash	125
Titanium	40
Calcium	5
Chlorine	30
Aluminum	20

RS Identifier	RS Description	Risk
	Cobalt: No. CAS 7440-48-4 EC No. 231-158-0	This device or one or more components of the device contains the substance in question, defined as CMR 1B, at a concentration above 0.1% by weight. Current scientific evidence and biocompatibility studies support that medical devices manufactured from cobalt alloys or stainless steel alloys do not cause an increased risk of cancer or adverse reproductive effects.

WARNINGS

1. The packages should be checked for tears or damage before the operation. If the sterile barrier is damaged, send the product back to Zimed Medical.
 2. Hip Replacement products are presented as sterile. It is not appropriate to re-sterilize the products. Please notify the manufacturer of the products whose sterility is impaired.
 3. Devices are for single use only. Do not reuse the used product.
 4. Please check the expiration date. In case of expired product, please contact the manufacturer.
 5. Surgical instruments used in conjunction with the products may have been stressed by previously applied pressures, which may have created defects that could cause the device to malfunction. Therefore, the instruments should be inspected for wear and damage before use. Damaged surgical instruments should never be used.
 6. The surgeon should apply the application methods related to the implant as given in the surgical technique. If the surgical technique is taken into consideration, the implant will be used safely and effectively.
 7. Zimed Medical products have not been evaluated for safety and compatibility in MR environment.
 8. Interchangeable components according to EN ISO 21534 Annex C;
 - CoCrMo Alloy (ISO 5832-4)/ Titanium Alloy (ISO 5832-3)
 - CoCrMo Alloy (ISO 5832-4)/ CoCrMo Alloy (ISO 5832-4)
 - Implant Steel (ISO 5832-1)/Titanium Alloy (ISO 5832-3)
 - Implant Steel (ISO 5832-1)/Implant Steel (ISO 5832-1)
- Components that are not suitable for use with each other;
- Implant Steel (ISO 5832-1)/ CoCrMo Alloy (ISO 5832-4)

EXPECTED CLINICAL BENEFITS

TOTAL HIP REPLACEMENT

Total hip arthroplasty is expected to improve patient mobility and reduce pain by replacing the damaged hip joint articulation in skeletally mature patients with evidence of sufficient intact bone to replace and support the components.

PARTIAL HIP PROSTHESIS

Partial hip replacements are expected to improve patient mobility and reduce pain in skeletally mature patients with evidence of a satisfactory native acetabulum, adequate bone setting and replacement of the damaged hip joint articulation supporting the femoral stem.

ZİMED MEDİKAL

FEMORAL STEMS USER GUIDE

Attention of the Specialist (EN)

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Three different retrospective studies include comprehensive analyses evaluating various prostheses used in hip arthroplasty in terms of both safety and performance criteria. These analyses addressed safety parameters, including revision rate, infection rate, and dislocation rate, as well as performance parameters, including the Harris Hip Score (HHS) and mortality rate. Each of these retrospective studies, conducted at Harran University Hospital and Tarsus State Hospital, showed that the prostheses used in hip arthroplasty achieved successful results in terms of both safety (revision, infection, dislocation rates) and performance (HHS and mortality rates), and were highly consistent with the literature and current state-of-the-art (SOTA) criteria. Articles in the literature investigating the use of Zimed Medical femoral stems were reviewed, and three retrospective case series methodologies were included in the study. In these studies, it was observed that the hip prostheses used achieved successful outcomes in terms of safety and performance criteria, were consistent with the literature, and complied with the criteria defined in SOTA.

Information on safety and clinical performance summaries (SSCP) can be accessed at www.zimed.com.tr and <https://ec.europa.eu/tools/eudamed/#/screen/search-eo>.

USAGE INFORMATION

Certain methods may change over time as more clinical experience is gained. Surgical techniques can be accessed at www.zimed.com.tr.

Before Operation

1. The following factors should be considered in patient selection: Patient weight, activity level and occupation, comorbidities, alcohol, smoking and substance use, etc. Product life and durability can be affected by these variables. An overweight patient may place excessive loads on the product, leading to product failure.
2. Selection of a prosthesis suitable for the type of fracture is extremely important. The specialist should select the appropriate prosthesis according to orthopedic practices before the operation. Incorrect selection or positioning of components may affect performance by causing leg length discrepancy, limping, loosening, collapse, femoral impingement and periarticular calcification.
3. The specialist should warn the patient that the prosthesis cannot replace normal healthy bone, that the prosthesis may be broken or damaged by certain movements or trauma, that the expected duration of use is limited and that it may need to be replaced in the future.
4. Take extreme care when handling and storing the implant. If the surface of the implant is bent, cut or scratched, the implant may significantly reduce its fatigue resistance, strength resistance and/or wear characteristics. This can translate into invisible internal stress, which can lead to product breakage.
5. Information on surgical techniques is available upon request. The specialist should be familiar with the surgical technique.
6. The surgical instruments used with the implants may have been subjected to stresses due to previous stresses, which may have created defects that could cause the device to fail. Therefore, surgical instruments should be inspected for wear and damage before use. Damaged surgical instruments should never be used.
7. Before surgery, the packages should be checked for tears or damage. If the sterile barrier is broken, send the product back to Zimed Medical.

Operation

1. The operation should be performed in accordance with orthopedic conventions and taking into account the shapes specified in the surgical technique.
2. Care should be taken with regard to complications and the position of the prosthetic components in the bone should be checked by fluoroscopy.
3. During surgical intervention, surgical instruments may break and crack. Surgical instruments that are overused or subjected to excessive force may break. Hand instruments should be checked for wear and fractures before surgery.

After Operation

1. The patient should be informed of the limitations of the reconstruction and the need to protect the prosthesis from bearing the full weight until adequate fixation and healing is achieved. The patient should be cautioned to limit his/her activities and to protect the replaced joint from unreasonable stresses and possible loosening, fracture and/or abrasion and to follow the doctor's instructions regarding treatment.
2. The patient should be warned that the prosthesis is not a substitute for normal healthy bone, that the prosthesis may be fractured or damaged by movement or trauma, that its expected lifetime is limited and that it may need to be replaced in the future.
3. Wear of the components may lead to increased loosening and damage to the bone.
4. Periodic, long-term monitoring is recommended to monitor the position and condition of the prosthetic components and the condition of the adjacent bone. Periodic postoperative x-rays are recommended for close comparison with early postoperative conditions to find evidence of long-term changes related to displacement, loosening, bending or cracking of components.

SIDE EFFECTS

1. Cardiovascular disorders, including venous thrombosis, pulmonary embolism or myocardial infarction may occur.
2. Damage to blood vessels or hematoma may occur.
3. Delayed wound healing; deep wound infection (early or late) may occur, which may require removal of the prosthesis. In rare cases, arthrodesis of the involved joint or amputation of the limb may be required.
4. Osteolysis (progressive bone loss). Osteolysis can be asymptomatic and therefore routine periodic radiographic examination is vital to prevent future serious complications.
5. Unintended shortening or lengthening of the limb may occur.
6. Fatigue fracture of prosthetic components may occur as a result of trauma, strenuous activity, misalignment, incomplete implant placement, time in service, loss of fixation, nonunion or excessive weight.
7. Trochanteric nonunion may occur due to inadequate reattachment and/or premature weight bearing, resulting in prosthesis collapse.
8. Trochanteric avulsion may occur as a result of excessive muscle tension, premature weight loss or inadvertent intraoperative weight loss.
9. Temporary or permanent nerve damage, peripheral neuropathies and subclinical nerve damage may occur as a possible consequence of surgical trauma resulting in pain or numbness in the affected limb.
10. Sensitization to metal or histologic or allergic reaction or adverse foreign material reaction may occur due to implantation of a foreign material.
11. There may be tissue reactions to implants, including macrophage and foreign body reaction around the implant.
12. Inadequate range of motion may occur due to improper selection or positioning of components, femoral compression and periarticular calcification.

REMOVAL OF THE IMPLANT

The length of time the implant remains in the patient may vary depending on the patient's daily movements. The product life for the implant is set at 15 years. Implants can be removed according to the decision of specialist physicians.

PRODUCT DISPOSAL

Unless otherwise specified, implants should be disposed of as medical devices in accordance with healthcare facility procedures.

MR SAFETY INFORMATION

Zimed Femoral Stems have not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating or unintended movement in the MR environment. MR testing on a person with this medical device could result in injury or malfunction of the device.

STERILIZATION

Zimed Medical Femoral Stems products are individually packaged and supplied sterile.

In case of any negative situation with the product, please contact the manufacturer. Manufacturer information detail is given below.

PRODUCER: ZİMED Medikal San. ve Tic. Ltd. Şti.
Aydınlar Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel : 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

ZİMED MEDİKAL

FEMORAL STEMS USER GUIDE

Attention of the Specialist (EN)

Document No: IFU.07.02/02
17.11.2025

Sterilization method is indicated on the label on the package. Minimum 25 kGy gamma ray is applied to all products sterilized with gamma ray. Sterile products remain sterile until the expiry date unless package integrity is compromised. Re-sterilization of the products is not appropriate. Please notify the manufacturer if there is visible damage to the sterile packaging or if the packaging has been opened accidentally.

SHELF LIFE

Shelf life is 10 years for gamma sterilized Femoral Stems and 5 years for those sterilized with ethylene oxide.

Please check the expiration date. In case of expired product, please contact the manufacturer.

STORAGE AND STORAGE CONDITIONS

Implants should be stored in a dry, clean environment, protected from direct sunlight, extreme heat and humidity.

PRODUCT SHIPMENT

Sterile Femoral Stems are packaged in accordance with the sterilization method. All products are shipped shrink wrapped in a box with user manual and label.

SERIOUS INCIDENTS

Serious events occurring in connection with the prosthesis must be reported to the manufacturer and the competent authority of the Member State where the user and/or patient is located.

LABEL INFO

Labels are prepared with reference to Medical Device Regulation (MDR) 2017/745 and ISO 15223. A table with explanations about the symbols on the label and box is given below. The country code according to the language of the user manual is written with reference to the ISO 639-1 standard.

(See Symbol Glossary)

SYMBOL GLOSSARY		
SYMBOL	STANDARD REFERENCE	EXPLANATORY TEXT
	-	Company Logo
	European Medical Device Directive 93/42/EEC (as amended by Directive 2007/47/EC), Article 17 European Medical Device Regulation 2017/745 Madde 20	"0476" Notified Body Number
	ISO 15223-1 (Article 5.7.7*)	Medical Device
	ISO 15223-1 (Article 5.1.5) ISO 7000 (Symbol No. 2492)	Lot No
	ISO 15223-1 (Article 5.1.6) ISO 7000 (Symbol No. 2493)	Product Reference/Catalog Number
	ISO 15223-1 (Article 5.7.10*)	UDI number (Unique Device Identification Number)

	ISO 15223-1 (Article 5.4.3) ISO 7000 (Symbol No. 1641) e-IFU ISO 15223-1 Ek A / A.15	Electronic User Manual
	ISO 15223-1 (Article 5.1.1) ISO 7000 (Symbol No. 3082)	Manufacturer Information
	ISO 15223-1 (Article 5.1.3) ISO 7000 (Symbol No. 2497)	Manufacture Date
	ISO 15223-1 (Article 5.1.4) ISO 7000 (Symbol No. 2607)	Expiration Date
	ISO 15223-1 (Article 5.2.12*) ISO 7000 (Symbol No. 3704)	Double Sterile Barrier System
	ISO 15223-1 (Article 5.2.4) ISO 7000 (Symbol No. 2502)	Sterilized using irradiation
	ISO 15223-1 (Article 5.2.4) ISO 7000 (Symbol No. 2501)	Sterilized using ethylene oxide
	21 CFR 801.15(c)(1)(i)F 21 CFR 801.109	Use as directed by prescription
	ISO 15223-1 (Article 5.4.2) ISO 7000 (Symbol No. 1051)	Implants are for single use only
	ISO 15223-1 (Article 5.2.6) ISO 7000 (Symbol No. 2608)	Do not re-sterilize
	ISO 15223-1 (Article 5.2.8) ISO 7000 (Symbol No. 2606)	Do not use damaged packaging
	ISO 15223-1 (Article 5.3.2) ISO 7000 (Symbol No. 0624)	Keep away from direct sunlight
	ISO 15223-1 (Article 5.3.4) ISO 7000 (Regulation No.0626)	Keep away from rain
	ISO 15223-1 (Article 5.3.1) ISO 7000 (Symbol No. 0621)	Kırılabilir dikkatli taşıyınız
	ISO 15223-1 (Article 5.4.4) ISO 7000 (Symbol No. 0434A)	Warning
	ISO 15223-1 (Article 5.4.10)	Contains hazardous substances

Note Bold-italicized text are special warnings that the user should pay special attention to.