www.zimed.com.tr

zimed®

**İSKELET PİN VE TELLERİ
KULLANIM KILAVUZU**
*INSTRUCTIONS FOR USE
FOR SKELETAL PINS*

TR

ENG

ES

IT

FR

RU

innovation to move

(TR) ÜRÜN TANIMI:

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(EN) PRODUCT DEFINITION:

Skeletal Pins and Wires are designed for the open or percutaneous fixation of fractures, bone reconstruction, as guide pins for placement or assisting to the fixation of other implants, for percutaneous implantation and as traction pins applicable to the skeletal system.

Steinmann Nails are nails with cylindrical bodies and compact structure designed for use in temporary traction for bone fixation. It has 3 length options as 200, 250, and 300 mm at Ø3 mm, and 1 option as 300mm at Ø3,5; while there are 2 length options as 150 mm and 200 mm at Ø4 mm, and 3 length options as 150, 200, and 250 mm at Ø 4,5 mm. At Ø5 mm, there are 5 length options in total between 150 mm and 350 mm with increments of 50 mm. At Ø6 mm, there are 4 length options as 200, 250, 300 and 350 mm. Steinmann nail is presented as non-sterile. It is made from implant steel complying with the ISO 5832-1 standard.

Kirschner wire has 4 different types including trocar tip, bayonet tip, diamond Kirschner wire and threaded Kirschner wire used in the fixation of bone fractures. The trocar-tipped Kirschner wire has 2 different length options as 150 mm and 300 mm at Ø 0,8 size, 4 different length options as 250, 275, 300, and 400 mm at Ø 1,0 mm, and 3 different length options as 250, 300, and 400 mm at Ø 1,2 mm. At Ø 1,4 mm, there is 150 mm length option, and at Ø 1,5 mm, there are 5 different length options as 150, 250, 275, 300, and 400 mm. At Ø 1,8 mm, there are 3 different length options as 275, 300, and 400 mm. At Ø 2 mm, there are 6 different length options in total as 200, 250, 275, 300, 400, and 440 mm. At Ø 2,2 mm, there is 300mm length option. At Ø2,5 mm, there are 4 different length options in total as 250, 275, 300, and 400 mm. At Ø 2,6 mm, there is one single length of 300 mm. At Ø 3,0 mm, there are 3 different length options as 275, 300, and 400 mm. At Ø 3,5mm, there are 2 different length options as 250 mm and 300 mm, and at Ø4,0 mm, there is one length of 300 mm. Bayonet-tip Kirschner wire has 4 different length options as 300 mm and 400 mm at Ø 1,8 ad Ø 2,0 mm. The threaded Kirschner wire comes with length options of Ø 0.5 X 250 mm, Ø 0.5 X 300 mm, Ø 0.8 X 250 mm, Ø 0.8 X 300 mm, 2 options of Ø 0.8 X 350 mm and 400 mm, 4 different options at Ø 1.0 as 250 mm, 300 mm, 350 mm, and 400 mm, one length option of Ø 1.2 X 300 mm, two length options of Ø 1.5 X 150 mm and 300 mm, and 4 length options of Ø 1.8 X 250 mm, 300 mm, 350 mm, and 400 mm. There are 6 options of Ø 2.0 X 200 mm, 250 mm, 280 mm, 300 mm, 350 mm, and 400 mm. There is one length option of Ø 2.2 X 300 mm. There are 4 length options of Ø 2.5 X 250 mm, 300 mm, 350 mm and 400 mm. There are 3 length options of Ø 3.0 X 250 mm, Ø 3.0 X 300 mm, 350 mm, and 400 mm. There are 4 length options of Ø 3.5 X 250 mm, 300 mm, 350 mm, and 400 mm. There are 4 different length options in total as Ø 4.0 X 250 mm, 300 mm, 350 mm, and 400 mm. Diamond-tip Kirschner wire has one length of Ø2x440 mm. Kirschner wire is used as a guide for plaques. At the same time, it can act to assist fixation for external fixators, or can be used alone in small fractures. Kirschner wires are made from implant steel complying with the ISO 5832-1 standard and put into the market in non-sterile form.

Cable is manufactured to assist to the fixation of other implants in bone reconstruction. It is indicated in peri-prostatic fractures, re-fixing the trochanter, prophylactic cerclage, stabilization of bone graft material, fixation of allograft support, stabilization of fractures in primary or revision arthroplasty process and traumatic fractures requiring both coating and cerclage. The cable system comes in different measures including Ø 1,5 and Ø 2,0 mm, in lengths of 500 and 700 mm. It is made from Ti6Al4V, ISO 5832-1 implant steel and ASTM F 1537 CoCr raw material complying with ASTM F 136 standard. Cable materials are offered as non-sterile.

Cerclage Wires: These are non-fragile and soft connectable products used to connect bones after fractures. Cross-sectional area is varied between 0.5 and 1,5 mm, and the length is 200 mm. Cerclage wires are made from implant steel complying with ISO 5832-1 standard. It is put into the market as non-sterile.

U-nail is used for the fixation of connective tissue, bones or tendons and fixation of anatomic components. It is made from TiGr3 raw material complying with ISO 5832-2 standard. There are 6 different length options as 6, 8, 9,5, 10, 11 and 12 mm. U-nail is presented in both sterile and non-sterile forms. Sterile products are sterilized with gamma or EO sterilization. U-nail has different models including plain (punch) and graded models. These models have 3 sizes as S, M, and L.

ZİMED MEDİKAL

INSTRUCTIONS FOR USE FOR SKELETAL PINS

Submitted to the Attention of Surgeons (EN)

Document No: IFU.06/00
14.06.2022

INDICATIONS:

Skeletal Pins and Wires are used for the following purpose:

- To act as a guide for other implants when stabilizing fractures and/or placement in the skeletal system,
- To be used as a supplement in external fixation,
- To act as auxiliary fixator for Skeletal Pins and Wires,
- For bones suitable for the treatment in osteotomy, arthrodesis, fracture fixation and external fixation,
- Implantation in accordance with the principles of AO / ASIF,
- To act as a guidewire for implants in osteosynthesis treatment.

Do not use the product rather than the indications stated here.

PATIENT POPULATION Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

CONTRA-INDICATIONS:

Active infections, manifest or occult infection or manifest local inflammation in or around the affected site, existing fracture fixation device implanted before, previous bone deformation, hypovolemia, hypothermia, coagulopathy, obesity, advanced osteoporosis, poor bone quality or amount, disorders in the medullar canal, and poor blood circulation are among the contra-indications.

WARNINGS:

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
3. Products of different manufacturers must not be used together.
4. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
5. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning, and must inform the patient.
6. These products are for single use only. Do not re-use the product.

CONDITIONS THAT MAY AFFECT THE INTENDED PERFORMANCE:

Preoperative:

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.
2. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.
3. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension, and can result in breaking of the product.
4. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.
5. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must

be examined for wear and tear and damages. Damaged surgical hand tools must never be used.

6. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

Operative

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.

2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.

3. During the implantation, fractures can occur in the application area because of trauma or overloading, particularly in the presence of poor bone structure.

4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.

2. Abrasion of components can lead to increase of loosening and damaging of the bone.

3. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending of cracking of components.

4. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.

5. Recommendations of surgeons in the postoperative period are of utmost importance.

6. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

ADVERSE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism and myocardial infarct can occur.

2. Damaging of blood vessels or hematomas can be seen.

3. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.

4. Unwanted lengthening or shortening of the extremity can occur.

5. Fatigue fractures of products can arise from trauma, exhausting physical activity, misalignment, defective implantation, service period, loss of fixation, nonunion or overloading.

6. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.

7. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydinlar Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

8. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.

9. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT: The period of the implant within the patient's body is minimum 6 weeks and can be extended to 1 year or more in difficult situations; afterwards, it must be removed based on the surgeon's decision.

STERILIZATION:

Implants Put in the Market as Sterile: Skeletal pins and wires are put in the market as sterile or non-sterile. Sterile products will remain so unless the integrity of the packaging is broken. Two types of sterilization methods, gamma irradiation and ethylene oxide, are applied. Minimum 25 kGy gamma irradiation dosage is applied to all the implants sterilized with gamma irradiation.

Please notify the manufacturer about the product with impaired sterility if put into the market as sterile.

Implants Put in the Market as Non-Sterile: Products put in the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in.

All the skeletal pins and wires and surgical hand tools must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

Autoclave Sterilization: Non-sterile products put in the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

AUTOCLAVE STERILIZATION	
Parameters	Autoclave with pre-vacuum
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying Time	15 minutes

Note: Exposure time means the period to be applied starting from the moment that the autoclave temperature reaches the exposure temperature.

The table below shows autoclave types for specific sterilization parameters.

PRODUCT TYPE	ÖNERİLEN STERİLİZASYON METODU
Skeletal Pins and Wires	High Temperature Prevac Steam (Autoclave with Pre-Vacuum) Flash Prevac Steam (Pre-Vakum Autoclave with Pre-Vacuum) If used, practice must be only according to ANSI/AAMI S79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)

Since the method of use will change for each device, please consult the written instructions of your producer for your sterilizer and loading configuration.

SHELF LIFE:

Shelf life for **Skeletal Pins and Wires** sterilized with gamma irradiation is 5 years; while the same for those sterilized using ethylene oxide is 3 years.

Please check the expiry date. Contact the manufacturer for products expired.

Non-sterile products do not have shelf lives.

DISPATCH OF PRODUCTS:

ZİMED MEDİKAL

INSTRUCTIONS FOR USE FOR SKELETAL PINS

Submitted to the Attention of Surgeons (EN)

Skeletal Pins and Wires presented as sterile or non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrank together with instructions for use and label.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity. The recommended storage conditions for sterile products are 15-35°C.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the manufacturer and used and/or authorized entity of the member state that the patient is a citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

(See: Descriptions of Symbols)

DESCRIPTION OF SYMBOLS	
	Logo of the Firm
	"1984" Notified Body Number
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI Number (Single Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Expiry Date
	Double Sterile Barrier System
	Sterilized with irradiation
	Sterilized with ethylene oxide

Document No: IFU.06/00
14.06.2022

	Use with prescription
	Implants are for single use
	Do not re-sterilize
	Do not use the damaged package
	Keep away from direct sunlight
	Keep away from precipitation
	Fragile: Handle with care
	Temperature limits
	Caution
	Non sterile

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydınlı Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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