www.zimed.com.tr

zimed®

MİNİ PLAK SİSTEMİ
KULLANIM KILAVUZU
MINI PLATE SYSTEM
INSTRUCTIONS FOR USE

TR

ENG

ES

IT

FR

RU

innovation to move

(TR) ÜRÜN TANIMI:

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ZİMED MEDİKAL

MINI PLATE SYSTEM INSTRUCTIONS FOR USE

Submitted to the attention of Surgeons

(EN)

Document No: IFU.01.03/00
14.06.2022

(EN) PRODUCT DEFINITION:

It includes plates, screws, and hand tools that will enable the application of these products to ensure the anatomical alignment of the broken bones, to keep them together, and to provide splinting with the ability to ensure the early mobility of the patient. The products are non-sterile and the implants in the set are for single use. See product catalog for reference numbers and names of products.

Mini Plate System is made of titanium alloy (ISO 5832-2, ISO 5832-3, EN ISO 14602).

Locked LCP Straight Plate: It consists of 2.0 mm, 2.4 mm and different sizes designed plates in different shapes compatible with the anatomical structure of the region, for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery. There are 4,5,6,7,8 hole options for Ø2.0 mm. It is used together with Ø2.0 mm locked and Ø2.0 mm cortex screw. There are 4, 5, 6, 7, 8, 9, 10 hole options for Ø2.4 mm. Ø2.4 mm locked and Ø2.4 mm cortical screw are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP T Plate Type I: It is designed 2.0 mm in size for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery, in different shapes compatible with the anatomical structure of the region, For Ø2.0 mm, there are 2+3 (head 2 holes), 2+5, 2+7 holes options. It is used with Ø2.0 mm locked and Ø2.0 mm cortex screw. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP T Plate Type II: It is designed 2.0 mm in size for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery, in different shapes compatible with the anatomical structure of the region. For Ø2.0 mm, there are 2+3 (head 2 holes), 2+5, 2+7 hole options. It is used with Ø2.0 mm locked and Ø2.0 mm cortex screw. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP Reconstruction Plate: It consists of 2.0 mm and 2.4 mm different sizes designed plates in different shapes compatible with the anatomical structure of the region, for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery. There are 6, 12 holes options for Ø2.0 mm. It is used with Ø2.0 mm locked and Ø2.0 mm cortex screw. There are 6, 12 holes options for Ø2.4 mm. Ø2.4 mm locked and Ø2.4 mm cortical screw are produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP Y Plate: It consists of 2.0 mm, 2.4 mm and different sizes designed plates in different shapes, compatible with the anatomical structure of the region, for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery. There are 2+3 (head 2 holes), 2+5, 2+7 holes options. Ø2.0 mm plates are used with Ø2.0mm locking and Ø2.0 mm cortex screws, and Ø2.4 mm plates are used with Ø2.4mm locking and Ø2.4 mm cortical screws. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP T Plate: It consists of 2.0 mm and 2.4 mm different sizes designed plates in different shapes compatible with the anatomical structure of the region, for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery. There are 3+3 (head 3 holes), 3+5, 3+7 holes options. Ø2.0 mm plates are used with Ø2.0mm locking and Ø2.0 mm cortex screws, and Ø2.4 mm plates are used with Ø2.4mm locking and Ø2.4 mm cortical screws. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP T Phalange Plate Type I: It is designed for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery, in different shapes compatible with the anatomical structure of the region, with a size of 2.4 mm. For Ø2.4 mm, there are 2+3 (head 2 holes), 2+5, 2+7 holes options. It is used with Ø2.4 mm locked and Ø2.4 mm cortex screw. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked LCP T Phalange Plate Type II: It is designed for fusion, osteotomy and fracture determination in hand and foot orthopedic surgery, in different shapes compatible with the anatomical structure of the region, with a size of 2.4 mm. For Ø2.4 mm, there are 2+3 (head 2 holes), 2+5, 2+7 holes options. It is used with Ø2.4 mm locked and Ø2.4 mm cortex screw. It is produced from ISO 5832-2 TiGr3 (ASTM F 67) material.

Locked Screw: There are Ø 2.0 mm and Ø 2.4 mm options. It is used for fixing mini plates. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 6-30 mm.

Cortex Screw: There are Ø 2.0 mm and Ø 2.4 mm options. It is used for fixing mini plates. It is produced from Ti6Al4V ELI (ASTM F 136) material. Various length options are available between 6-30 mm.

Do not use the product rather than the indications stated here.

PATIENT POPULATION: Our products are designed in several sizes with the purpose of ensuring the geometric conformity between the fixator and the bone in patients with different anthropometric measurements. It is intended for universal use. General principles are applied in the selection of patients and implants. Selection of the correct implant is very important. The proper type and measure must be determined by taking anatomic and biomechanical factors including age, level of mobility, weight, conditions of bones and muscles of the patient, and if he/she had previous operations or not etc. into consideration.

CONTRAINDICATIONS: Sepsis, malignant primary or metastatic tumors, material sensitivity, damaged vasculature, pertrochanteric fractures or trochanteric fractures with subtrochanteric enlargement (31-A3), overweight patients, patients with insufficient skin tissue, bone or neurovascular condition, damaged tendon, patients with high activity levels, patients with diabetes, patients with osteoporosis, patients with blood circulation problems, patients with material sensitivity, psychologically or physiologically incompetent patients, elderly patients, active and latent infection

1. Zimed Medikal implants may only be used by authorized specialists and be implanted under sterile operation room conditions.
2. In case of suspected hypersensitivity in the patient, tests must be carried out by all means before the implantation.
3. Products from different manufacturers must not be used together.
4. Unless otherwise stated, implants must be disposed of as medical devices according to the procedures of the healthcare facility.
5. Patients can be hypersensitive or allergic against implant materials, and metal ions in particular. The type of material is shown on the box label for each implant. The surgeon must evaluate the possible risks related to implant materials as the part of the preoperative planning and must inform the patient.
6. These products are for single use only. Do not re-use the product.
7. The products are anatomically designed. Plate bender should not be used for anatomical products.

SITUATIONS THAT MAY AFFECT INTENDED PERFORMANCE:

Preoperative:

1. The following factors, which can increase the risk of failure and can be critically important for the final outcome of the procedure must be taken into consideration in patient selection: weight, activity level and occupation of the patient, smoking or alcohol use etc. Product life and durability can be affected by these variables. An overweight patient can load excessive pressure on the product resulting it to fail.
2. Surgeon must also inform the patient about any other risks she/he can deem required to be acknowledged to the patient.
3. Please pay utmost attention to the handling and storage of the implant. Twisting, cutting or cracking of the implant surface can significantly reduce the fatigue resistance, strength and/or erosion characteristics of the implant. These can transform into invisible internal tension and can result in breaking of the product.
4. These can transform into invisible internal tension and can result in breaking of the product.
5. Information on procedural techniques is provided upon request. The surgeon must be familiar with procedural techniques.
6. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be

examined for wear and tear and damages. Damaged surgical hand tools must never be used. The surgical hand tools used with the products can have defects arising from pressure in previous uses that can lead to failure of the product. Therefore, surgical hand tools must be examined for wear and tear and damages. Damaged surgical hand tools must never be used. Packages must be checked for tears or damages in the preoperative period. Please return the product to Zimed Medikal if the sterile barrier has been broken.

7.

Operation

1. Known surgical principles and for the selection of patients, general principles must be followed. Selection of the right implant has utmost importance. Any past or future operations that the patient had undergone/may undergo and anatomical and biomechanical factors including age, activity level, weight, and bone and muscle conditions of the patient must be considered when selecting the type and measures that will fit the patient. In general, component with the largest section area that will provide adequate support will be preferred.
2. The surgeon must use the application methods relevant for the implant as indicated in the surgical techniques section. The implant applied according to the surgical technique will be used effectively and safely.
3. During the implantation, fractures can occur in the application area because of trauma or overloading, particularly in the presence of poor bone structure.
4. During the surgical intervention, surgical hand tools can break or crack. Much-used surgical hand tools or those used by applying excessive pressure can break. Surgical hand tools must be checked for abrasion or breaks before the operation.

Postoperative

1. The patient must be informed about the limitations of reconstruction and the requirement of avoiding to load the product with full weight of the patient until adequate fixation and full recovery. The patient must be warned about the necessity of limiting his/her activities and protecting the replaced joint from unreasonable pressure and potential loosening, fractures and/or wear and tear and to follow the doctor's instructions concerning the treatment.
2. The patient must be warned that the implant cannot replace a normal and healthy bone, the implant can break or be damaged as a result of some movements and traumas, its expected life is limited, and replacement can be necessary in the future.
3. Abrasion of the components can result in the increase of loosening and damaging of the bone.
4. Long-term periodic follow-up is recommended to observe the positions of prosthetic components and condition of the adjacent bone. Postoperative periodic X-rays are recommended for close comparison with the purpose of finding evidences of long-term changes related to displacement, loosening, bending of cracking of components.
5. Zimed Medikal products have not been evaluated for safety and compliance in MR environment.
6. Recommendations of surgeons in the postoperative period are of utmost importance.
7. Postoperative treatment must be arranged so as to ensure gradual recovery of muscular strength in the application area.

SIDE EFFECTS:

1. Cardiovascular disorders including venous thrombosis, pulmonary embolism or myocardial infarcts can occur.
2. Damaging of blood vessels or hematomas can be seen.
3. Osteolysis (progressive bone lysis). Osteolysis can be asymptomatic, and therefore, periodic radiographic examinations are vitally important to prevent any possible severe complications in the future.
4. Unwanted lengthening or shortening of the extremity can occur.
5. Fatigue or breaks of products can happen as a result of trauma, exhausting activity, wrong alignment, defective placement of implants, period of service, loss of fixation, lack of union or excessive weight

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydinlar Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 **Fax:** 0 342 238 44 11 **Web:** www.zimed.com.tr **E-mail:** info@zimed.com.tr

ZİMED MEDİKAL

MINI PLATE SYSTEM INSTRUCTIONS FOR USE

Submitted to the attention of Surgeons
(EN)

Document No: IFU.01.03/00
14.06.2022

6. As a potential result of surgical trauma resulting in pain or numbness in the affected extremity, temporary or permanent nerve damage, peripheral neuropathies or subclinical nerve damages can be seen.
7. Sensitivity against the metal or histological or allergic reactions or adverse foreign material reaction can occur as a result of implantation of a foreign material.
8. Trochanteric non-union may occur due to inadequate re-attachment and/or premature weight bearing.
9. Trochanteric avulsion may occur as a result of excessive muscle tension, premature weight loss, or accidental intraoperative attenuation.
10. Macrophages against the implant or tissue reactions containing foreign material reactions can be seen around the implant.
11. Embolism can occur in patients who have extreme sedentary lifestyle.

REMOVAL OF THE IMPLANT:

Implants Put on the Market as Non-Sterile:

Products put on the market without sterilization must be sterilized before using. Such products must be sterilized in the hospitals they are used in. all Mini Plate System products and surgical instruments used with them must be accepted as non-sterile unless clearly labeled as sterile. The method of sterilization required for the products must be determined before the operation.

Autoclave Sterilization: Non-sterile products put on the market without sterilization must be sterilized before using. The recommended sterilization for these products is autoclave sterilization. Parameters of Autoclave Sterilization are given in the table below.

AUTOCLAVE STERILIZATION	
Parameters	Pre-vacuumed autoclave
Exposure Temperature	134°C
Exposure Time	5 minutes
Drying time	15 minutes

Note: Exposure time is the period to be applied starting from the moment that the autoclave reaches the exposure temperature.

The following table shows the autoclave types as compared according to certain sterilization parameters.

TYPE OF PRODUCT	RECOMMENDED STERILIZATION METHOD
Mini Plate System	High Temperature Prevac Steam (Pre-Vacuum Autoclave) Flash Prevac Steam (Pre-Vakum Otoklav): If used, practice must be only according to ANSI/AAMI ST79 provisions. (2010- Detailed guidelines for steam sterilization and sterilization assurance in healthcare facilities)

Since the method of use will change for each device, please consult the written instructions of your manufacturer for your sterilizer and loading configuration.

SHELF LIFE:

There is no shelf life for non-sterile products.

DISPATCH OF THE PRODUCTS:

Mini Plate System presented as non-sterile are packaged in accordance with the sterilization method. All the products are dispatched within the box as shrunk together with instructions for use and label.

STORAGE AND MAINTENANCE CONDITIONS:

Implants must be maintained in a dry and clean environment and protected from direct sunlight, excessive heat and humidity.

SEVERE EVENTS: Any severe events arising in relation with the implant must be notified to the producer and used and/or authorized entity of the member state that the patient is a

citizen of.

LABEL INFORMATION: Labels are prepared with reference to the Medical Devices Directive 93/42/EC and ISO 15223. The table including the explanations of symbols placed on the label and box is below. The instructions for use have been prepared with reference to the ISO 639-1 standard in the languages of the Community.

(See: Description of Symbol)

SYMBOL DESCRIPTORS	
	Logo of the Firm
	Number of the Notified Body "1984"
	Medical Device
	Lot Number
	Product Reference/Catalog Number
	UDI number (Single-Device Identification Number)
	Electronic Instructions for Use
	Details of the Manufacturer
	Date of Production
	Used only with prescription.
	Implants are for single use
	Do not use the damaged package
	Keep away from direct sunlight
	Keep away from precipitation
	Fragile: Handle with care

	Temperature limits
	Caution
	Non sterile

Note: Expressions written in bold italics are special warnings that the user must pay meticulous attention.

Please contact the manufacturer in case of any adverse event associated with the product. Details of the manufacturer are given below.

MANUFACTURER: Zimed Medikal San. ve Tic. Ltd. Şti.
Aydinlar Mah. 03070 Cad. No:4 Şehitkâmil 27580 Gaziantep-Türkiye
Tel: 0 342 238 43 44 Fax: 0 342 238 44 11 Web: www.zimed.com.tr E-mail: info@zimed.com.tr

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